

Technical Data

SLIMYONIK® AIR BODYSTYLER	
POWER SUPPLY	
Input voltage	230 V
Input frequency	50 Hz
Maximum power consumption	180 W
Fuse	2 A, T / 250 V
Tripping characteristic of the fuse	T (slow-blow), configuration: 5 × 20 mm
SAFETY	
Safety class	I
Degree of protection	IP20
MASSAGE PRESSURE RANGE	
Working pressure	2.67 kPa – 14.0 kPa
Minimum pressure	27 mbar
Maximum pressure	140 mbar
MISCELLANEOUS	
Dimensions	40 × 40 × 22 cm
Weight	9 kg

SLIMYONIK® Trousers	
Function	24-chamber compression trousers for lymph drainage 12 chambers on each side
Colour	Burgundy
Weight including tube	5 kg
Dimensions [cm*]	Maximum waist measurement 130 Length 152
Cleaning	Damp cloth; do not use any solvents or cleaning agents
Sanitary trousers	Not included in delivery – see Section 2.7, Accessories
Equipment	3 leg widths, adjustable using zip. Back of legs: 3 pressure compensation folds. Angled tube connector flap with hook-and-loop fastener.
Material pressure range	0–140 mbar
Pressure for cosmetic applications	MAX 80 mbar
Pressure for freely selectable algorithms	MAX 140 mbar
Expander	Not included in delivery – see Section 2.7, Accessories

* All cm measurements may have a tolerance of approx. 2 cm



SLIMYONIK® Jacket Not included in delivery	
Function	24-chamber compression jacket for lymph drainage 12 chambers on each side
Colour	Burgundy
Weight including tube	5 kg
Dimensions [cm*]	Sleeve length: 95, sleeve width: 16 at the wrist, 24 at the shoulder
Cleaning	Damp cloth; do not use any solvents or cleaning agents
Sanitary jacket	Not included in delivery – see Section 2.7, Accessories
Equipment	Front: 5 adjustable straps (17 cm* variation), 2 zips. Sleeves with zip: can be completely opened for cleaning. Back with zip: separates jacket into 2 halves. Tube cover flap with zip.
Material pressure range	0–140 mbar
Pressure for cosmetic applications and freely selectable algorithms	MAX 50 mbar
Recommended starting pressure	MAX 30 mbar
Expander	Not included in delivery – see Section 2.7, Accessories

* All cm measurements may have a tolerance of approx. 2 cm



2.7 Accessories	Description (see current price list)
Sanitary trousers	For use under compression trousers; use once for each treatment
Sanitary jacket	For use under compression jacket; use once for each treatment
Nasal cannula	With single-use antibacterial filter
Expander set for trousers	Used to extend the leg width by approx. 10 cm
Expander set for jacket	Used to extend the arm width and waist by approx. 6 cm

Manufacturer

beauty lumis GmbH
Bruckmannring38
85674 Oberschleissheim
Germany
Phone: +49 89 520 33670
Email: info@beautylumis.com
